

SIZE

Height:	87.6 cm (34.5 in.)
Width:	67.3 cm (26.5 in.)
Depth:	58.5 cm (23.0 in.)

BATH CAPACITY

With 8 vessels:	19 L (5 gal)
-----------------	--------------

SPINDLE SPEED

Range:	25 to 250 rpm
Accuracy:	±1 rpm
Displayed Resolution:	0.1 rpm

WEIGHT

Main unit, dry:	63.5 kg (140 lb)
Main unit, bath and vessels filled with water:	90.8 kg (200 lb)

TEMPERATURE

Programmable temperature range:	25.0 to 55.0 °C
Accuracy:	±0.5 °C
Control:	±0.1 °C

DIGITAL TEMPERATURE PROBES (ALL PROBES)

Range:	10 to 60 °C
Accuracy:	±0.1 °C from 30-55 °C
Displayed Resolution:	0.01 °C
Minimum Submersion Depth:	25 mm
Response Time:	<45 s (normally 30 s)

POWER

Input:	100-230 V ± 10%, 50-60 Hz, 1.5 A max
Output:	24 VDC, 3.75 A, 90 W



VISION G2 HEATER

Weight:	4.1 kg (9 lb)
Height:	22.9 cm (9.0 in.)
Width:	24.1 cm (9.5 in.)
Depth:	16.5 cm (6.5 in.)
Power Input:	100-230 V ± 10%, 50-60 Hz
Power:	700-1200 W

EXTERNAL POWER SUPPLY

Power Input:	100-230 V ± 10%, 50-60 Hz, 1.5 A max
Output:	24 VDC, 3.75 A, 90 W