

**SIZE**

<b>Height:</b>	67.3 cm (26.5 in.)
<b>Width:</b>	39.4 cm (15.5 in.)
<b>Depth:</b>	58.4 cm (23.0 in.)

**BATH CAPACITY**

<b>With 6 vessels:</b>	13.3 L (3.5 gal)
------------------------	------------------

**SPINDLE SPEED**

<b>Range:</b>	25 to 250 rpm
<b>Accuracy:</b>	±1 rpm
<b>Displayed Resolution:</b>	0.1 rpm

**WEIGHT**

<b>Main unit, dry:</b>	29.5 kg (65 lb)
<b>Main unit, bath and vessels filled with water:</b>	49 kg (108 lb)

**TEMPERATURE**

<b>Programmable temperature range:</b>	25.0 to 55.0 °C
<b>Accuracy:</b>	±0.5 °C
<b>Control:</b>	±0.1 °C

**DIGITAL TEMPERATURE PROBES (ALL PROBES)**

<b>Range:</b>	10 to 60 °C
<b>Accuracy:</b>	±0.1 °C from 30-55 °C
<b>Displayed Resolution:</b>	0.01 °C
<b>Minimum Submersion Depth:</b>	25 mm
<b>Response Time:</b>	<45 s (normally 30 s)

**POWER**

<b>Input:</b>	24 VDC, 3.75 A, 90 W
---------------	----------------------

Only use external power supply (listed to the right).



**VISION G2 HEATER**

<b>Weight:</b>	4.1 kg (9 lb)
<b>Height:</b>	22.9 cm (9.0 in.)
<b>Width:</b>	24.1 cm (9.5 in.)
<b>Depth:</b>	16.5 cm (6.5 in.)
<b>Power Input:</b>	100-230 V ± 10%, 50-60 Hz
<b>Power:</b>	700-1200 W

**EXTERNAL POWER SUPPLY**

<b>Power Input:</b>	100-230 V ± 10%, 50-60 Hz, 1.5 A max
<b>Output:</b>	24 VDC, 3.75 A, 90 W