

VISION[®] AUTOPLUS[™]

WEIGHT (DRY)

21.3 kg (47.0 lb)

SIZE

Height: 46.8 cm (18.4 in.)

Width: 38.1 cm (15.0 in.)

Depth: 66.0 cm (26.0 in.)

SAMPLER VOLUME ACCURACY

0.5 to 5.0 mL: ±0.05

5.1 to 10.0 mL: ±0.10

10.1 to 15.0 mL: ±0.15

POWER

Input: 24 VDC, 3.75 A, 90 W

Only use external power supply (listed below).

SAMPLING

**with an optimized method*

Minimum sample interval: 1 min* (5 min is normal)

Maximum sample interval: 999 h, 59 min (999:59)

Sample time points: 50 max

Maximum sample volume: 90% of vial capacity

COLLECTION

13 x 100 rack: 18 x 8 test tubes, 144 total

16 x 100 rack: 18 x 8 test tubes, 144 total

12 x 32 rack: 21 x 8 vials, 168 total

Waters UPLC[®] rack (12 x 32): 3 trays 6 x 8 vials, 144 total

Septum: Pre-cut type only



EXTERNAL POWER SUPPLY

Power Input: 100-240 V~, 50-60 Hz, 1.5 A max

Output: 24 VDC, 3.75 A, 90 W

VISION[®] AUTOFILL[™]**WEIGHT (DRY)**

14.5 kg (32.0 lb)

SIZE**Height:** 40.1 cm (15.8 in.)**Width:** 28.0 cm (11.0 in.)**Depth:** 47.3 cm (18.6 in.)**POWER****Input:** 24 VDC, 3.75 A, 90 W

Only use external power supply (see reverse).

