

SAMPLE FILTRATION

PATENT
PENDING

VISION[®] AUTOFILTER[™] CHANGER

The Vision[®] AutoFilter[™] Changer

automatically introduces 25 mm syringe filters into the fluidic path of your dissolution samples. Designed as a seamless addition to Hanson's Vision[®] family of quality test instruments, the Vision AutoFilter Changer brings a new approach to automated filtration by giving you the power to determine the life of your filters.

Rather than discard filters after one-time use, the Vision AutoFilter Changer enables you to automatically introduce filters to the fluidic path, remove them if desired, return them back into position, or completely discard them, all of which is programmable. There are almost no limitations to the Vision AutoFilter Changer — you decide what works best for your method.



VISION[®] AUTOFILTER[™] CHANGER FEATURES

- Filters 6, 7, or 8 samples simultaneously
- Automatic insertion and removal of filters from fluidic path
- Compatible with a variety of 25 mm syringe filters
- Supports automation-compatible filters
- Comes standard with both 18-filter and 25-filter stack tubes
- Customizable filter capacity (the ceiling is the limit!)
- Pull-through or push-through system configuration capability
- User can create, manage, and save up to 100 protocols
- Upload and download programs with flash drive
- User-defined alarms for sampling, maintenance, and calibration
- Ability to temporarily remove filters from the fluidic path for complete system rinse between samples
- Full color touchscreen with intuitive menus and programming
- Electronics and firmware include 32-bit digital technology
- Print to both serial and network (PostScript) printers
- Front-facing filter positioning for easy visual verification of filters
- Catch bin drawer conveniently stores discarded filters
- Elegant, ergonomic design with workhorse performance
- 21 CFR Part 11 security compliant
- USP / US FDA / ASTM / JP / EP / CE / CSA / WEEE / RoHS compliant



9810 Variel Avenue
Chatsworth, CA 91311 USA

Toll Free: 800.821.8165
Phone: 818.882.7266

Fax: 818.882.9470
www.hansonresearch.com

VISION® AUTOFILTER™ CHANGER



PUSH-THROUGH CONFIGURATION (Recommended)



PULL-THROUGH CONFIGURATION



ORDERING INFO

For more information, please contact your local Hanson representative or email us at:

sales@hansonresearch.com

43-101-008 Vision AutoFilter Changer, 115/230 V



9810 Variel Avenue
Chatsworth, CA 91311 USA

Toll Free: 800.821.8165
Phone: 818.882.7266

Fax: 818.882.9470
www.hansonresearch.com